

XYL

Xylazine Test (Urine) Dip Card

A rapid test for the qualitative detection of Xylazine in human urine.
For Forensic use only.

INTENDED USE

The **Xylazine Test (Urine) Dip Card** is a lateral flow chromatographic immunoassay for the detection of Xylazine in human urine.

Test	Calibrator	Cut-off
Xylazine (XYL)	Xylazine	50 ng/ml

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

SUMMARY

Xylazine (colloquially known as tranq/tranq dope) is a non-opioid tranquilizer used as a sedative, analgesic, and muscle relaxant in animals such as horses and cattle^{1,3}. In humans, it could cause central nervous system depression, respiratory depression, bradycardia, hypotension, and even death¹. Most of the non-fatal cases required medical intervention. In recent years, xylazine has emerged as an adulterant in recreational drugs such as heroin². Chronic use of xylazine is reported to be associated with physical deterioration and skin ulceration. Combining xylazine with other drugs that cause central nervous system depression compounds the sedative effects and can increase the risk of overdose and death².

PRINCIPLE

The **Xylazine Test (Urine) Dip Card** is a rapid chromatographic immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Xylazine, if present in the urine specimen below 50 ng/ml, will not saturate the binding sites of the antibody coated particles in the test dip card. The antibody coated particles will then be captured by immobilized Xylazine conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Xylazine level exceeds 50 ng/ml because it will saturate all the binding sites of anti-Xylazine antibodies.

A drug-positive urine specimen will not generate a colored line in the test line region, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test dip card contains mouse monoclonal anti-Xylazine antibody-coupled particles and Xylazine-protein conjugate. A goat antibody is employed in the control line system.

PRECAUTIONS

- For Forensic use only. Do not use after the expiration date.
- The test dip card should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in

the same manner as an infectious agent.

- The used test dip card should be discarded according to federal, state and local regulations.

STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30°C). The test dip card is stable through the expiration date printed on the sealed pouch. The test dip card must remain in the sealed pouch until use. **DO NOT FREEZE**. Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible particles should be centrifuged, filtered, or allowed to settle to obtain clear specimen for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For long-term storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

MATERIALS

Materials Provided

- Test dip card
- Package insert

Materials Required but Not Provided


- Specimen collection container
- Timer

DIRECTIONS FOR USE

Allow the test device, and urine specimen to come to room temperature [15-30°C (59-86°F)] prior to testing.

- Remove the test device from the foil pouch.
- Remove the cap from the test device. Label the device with patient or control identifications.
- Immerse the absorbent tip into the urine sample for 10-15 seconds. Urine sample should not touch the plastic device.
- Replace the cap over the absorbent tip and lay the device flatly on a non-absorbent clean surface.
- Read results at 5 minutes.

DO NOT INTERPRET RESULT AFTER 10 MINUTES.



INTERPRETATION OF RESULTS

(Please refer to the illustration above)

NEGATIVE: * **Two lines appear.** One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Xylazine concentration is below the detectable level (50 ng/mL).

***NOTE:** The shade of red in the test line region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Xylazine concentration exceeds the detectable level (50 ng/mL).

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test dip card. If the problem persists, discontinue using the lot immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory testing practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The **Xylazine Test (Urine) Dip Card** provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.^{4,5}
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- Test does not distinguish between drugs of abuse and certain medications.

PERFORMANCE CHARACTERISTICS

Reproducibility

Reproducibility studies were carried out using commercially available stork solutions of the drug analytes listed. Dilutions were made from the stork solution of each drug to the concentrations specified in the following tables. The results are listed in the following tables.

Xylazine conc.(ng/mL)	Total number of Determinations	Result	Precision
No Drug present	40	40 negative	>99%
25	40	40 negative	>99%
75	40	40 positive	>99%
100	40	40 positive	>99%

Analytical Sensitivity

A drug-free urine pool was spiked with drugs to the concentrations at $\pm 50\%$ cut-off and $\pm 25\%$ cut-off. The results are summarized below.

XYL Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
25	-50%	30	30	0
37.5	-25%	30	30	0
50	Cut-off	30	12	18
62.5	+25%	30	0	30
75	+50%	30	0	30

Analytical Specificity

The following table lists the concentration of compounds (ng/mL) that were detected positive in urine by The **Xylazine Test (Urine) Dip Card** at a read time of 5 minutes.

Drug	Concentration (ng/ml)
Xylazine	50
4-Hydroxy xylazine	50
±MBDB HCL	25,000
Clonidine	25,000
Butylone Hcl	780
N-Acetylprocainamide	70,000

Effect of Urinary Specific Gravity

Twelve (12) urine samples of normal, high, and low specific gravity ranges 1.000 to 1.035 were spiked with drugs at 50% below and 50% above cut-off levels respectively. The **Xylazine Test (Urine) Dip Card** was tested in duplicate using ten drug-free urine and spiked urine samples. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to pH ranges of 4.00 to 9.00 in 1 pH unit increment and spiked with the target drug at 50% below and 50% above cut-off levels. The spiked, pH-adjusted urine was tested with the **Xylazine Test (Urine) Dip Card**. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Xylazine positive urine. The following compounds show no cross-reactivity when tested with the **Xylazine Test (Urine) Dip Card** at a concentration of 100 µg/mL.

Non-Cross-Reacting Compounds

Citalopram HBr	Dexamethasone acetate	morphine
L-Thyroxine	Dextromethorphan	Mosapride Citrate
y-aminobutyric acid	Diazepam	Nalidixic acid
(±)-Octopamine hydrochloride	Diclofenac sodium	Nalorphine
(-)-11-nor-9-carboxy-delta	Dicyclomine	Naloxone hydrochloride
(-)-Cotinine	Diethylstilbestrol	Naltrexone Hydrochloride
(±)-Epinephrine hydrochloride	Diffunisal	Naproxen
(±)-Methamphetamine	Dihydrocodeine HCl	Narcotine
(±)-N'-Nitrosornicotine	Diltiazem	N-Benzylisopropylamine
(±)-Nornicotine	Dimenhydrinate	Nicotinamide
(1R,2S)-(-)-Ephedrine HCl	Diphenhydramine HCl	Nicotinic acid
(R,S)-Norcotinine	Dirithromycin	Nifedipine
(S)-(-)-Nicotine-d3	Doxypiramide	Nimodipine
(S)-N-Nitrosoanabasine (NAB)	Domperidone	Nitazene
(S)-N-Nitrosoanatabine	Dopamine Hydrochlorid	Norbuprenorphine
(t)-methadane	Doxylamine succinate	Nordoxepin hydrochloride
±MDEA	Droperidol	Norethindrone
1-Adamantanamine hydrochloride	Duloxetine	Norflloxacin
2-Phenylethylamine hydrochloride	Ecgonine	Normorphine
3,4-Methylenedioxypropavolone HCl	Ecgonine methyl Esten	Norsertaline -13C6 HCl
4CN-CUMYL-BUTINACA	EMDP HCl	Nortrityline Hcl

4F-ABUTINACA	Enalapril Maleate	Noscapine	Cannabidiol	JWH-073 4-Hydroxybutyl	Ranitidine
4F-MDMB-BUTICA	Enoxacin	Ofloxacin	Captopril	JWH-122 5- Hydroxypentyl meabolite	Rifampicin
5,5-Diphenylhydantoin	Erythromycin	O-Hydroxyhippuric acid	Carbamazepine	JWH-210 5-Hydroxypentyl meabolite	Risperidone
5F-ADBICA	Esomeprazole Magnesium	Olanzapine	Carisoprodol	JWH-250 4-Hydroxypentyl	S(-)Cathinone
5F-EMB-PICA	EStroline	Omeprazole	Cefaclor	JWH-250 5-Hydroxypentyl metabolite	S(+)-methamphetamine
5F-MDMB-PICA	Ethopropazine hydrochloride	Ondansetro	Cefadroxil	Ketamine	Salbutamol
5-Hydroxytryptamine	Ethylmorphine	Oxalic acid	Cefalexin	Ketoconazole	Salicylic
6-MAM	Ethylone	Oxazepam	Cefixime	Ketoprofen	S-Doxylamine
7-beta-Estradiol	Ethyl-p-aminobenzoate	Oxolinic acid	Cefoperazone sodium	Kynurenic	Serotonin
8-Chloro Caffeine	Fenofibrate	Oxycodone	Cefuroxime Axetil	Lactose	Sertraline HCl
AB-PINACA	Fenoprofen	Oxymorphone	Cephalexin	Lamotrigine	Sildenafil
4-Hydroxypentyl metabolite	Fexofenadine Hydrochloride	OXYPHENBUTAZONE	Cephadrine	Lansoprazole	Simvastatin
Acetaminophen	Fluconazole	Paclitaxel	Chloramphenicol	Levetiracetam	Sodium Valproate
Acetylsalicylic	Fluoxetine hydrochloride	p-Aminobenzoic	Chloroquine diphosphate salt	Levofloxacin Hydrochloride	Spice Cannabinoid Mix
Acyclovir	folic acid	Paroxetine	Chlorothiazide	Levonorgestrel	Spice Cannabinoid Mix 2
ADB-4en-PINACA	FUB-144	Penfluridol	Chlorpromazine hydrochloride	Levorphanol tartrate	Spice Cannabinoid Mix 3
Albumin	Furosemide	Penicillin	Cholesterol	Levothyroxine Sodium	Spironolactone
Amikacin	Gabapentin	Penicillin G Sodium Salt	Cimetidine	Lignocaine	Sulfamonomethoxine
Amiloride	Gatifloxacin	Pentachlorophenol	Citalopram hydrobromide	Lisinopril	Tadalafil
Amiodarone Hydrochlorid	Gemfibrozil	Pentobarbital	Citicoline Sodium	Lodeine	Telmisartan
Amitriptyline	Gentamicin	Perphenazine	Clarithromycin	Lomefloxacin Hydrochloride	Thebaine
Amlodipine besylate	Gentamicin phenethylamine	phenethylamine	Cloethocarb	Loperamide hydrochloride	Theophylline
Amoxicillin	Glibenclamide	Phenobarbital	Clomipramine	Loratadine	Thiamine
Ampicillin	Gliclazide	Phentermine	Clomipramine	Lorazepam gluronide	Tolazamide
APINACA (AKB-48)	Glipizide	Phenytol Sodium	Clopidogrel Bisulfate	Maprotiline HCl	Topiramate
5-Hydroxypentyl metabolite	Glucose	p-Hydroxymethamphetamine	Clopidogrel Hydrogen Sulfate	MDMA	Tramadol hcl
Apomorphine hydrochloride	Hemoglobin	Pioglitazone Hydrochlorid	clorprenaline	MDMB-4en-PINACA	Trazodone HCL
a-Pyrrolidonorale phenore hcl	Herion	Piracetam	Clozapine	Meperidine	Triamterene
Aripiprazole	Hydralazine	Prazosin hydrochloride	Cocaethylene	Metformin	Triazolam
ascorbic acid	Hydrochlorothiazide	Prednisone	Cocaethylene	Methaqualone	Trimethobenzamide Hydrochloride
Atenolol	Hydrocodone	Procainamide hydrochloride	Cocaine	Methoearbm01	Trimethoprim
Atomoxetine hydrochloride	Hydrocortisone	Procaine hydrochloride	Codeine	methoxyphenamine	Tryptamine
Atorvastatin	Ibuprofen	Promazine	Cortisone	Methylmorphinan	Tyramine hydrochloride
Atrazine	Iclilin	Promethazine hydrochloride	Creatinine	Metoclopramide hydrochloride	Ursodeoxycholic Acid
Atropine	Impipramine Hydrochloride	Propylthiouracil	d l 3-4-MDA	Metoprolol	Valproic Acid
Atroscine	Indapamide	Pseudoephedrine Hydrochloride	D L-Tryptophan	Metronidazole	Venlafaxine hydrochlorid
Azithromycin	Iproniazid	Pseudoephedrine(r)	D/L-Tyrosine	Midazolam	Verapamil hydrochloride
Baclofen	Isoproterenol hydrochloride	Pseudoephedrine(s)	Deoxycorticosterone	Mifepristone	Vitamin C
Benzocaine	Isoproterenol- (+/-)	Quetiapine fumarate	Desalkyl futazepam	Minocycline	Zomepirac
Benzoic acid	Isosorbide Dinitrate	quinidine	Desipramine HCl	Zonisamide	α-Hydroxyhippuric acid
Benzoyllecgonine	Isoxsuprine	Quinine	Desloratadine	Mitragynine	β-Estradiol
Berberine	JEH-018 5-Hydroxypentyl	Quinine Monohydrochloride Dihydrate	Desloratadine Citrate	MMB-CHMICA	R(-)Amphetamine
Bilirubin	JEH-018	R(-)-Methamphetamine (levo-Methamphetamine)	Disodium deta-9-THC	Montelukast sodium	Clobazam
Buprenorphine	N-(4-hydroxypentyl) metabolite	R(-)-Phenylephrine HC	6β-Naltrexol	Alprazolam	Hydrocodone
Bupropion hydrochloride	JWH-019 5-Hydroxyhexyl	R(+)-Methcathinone HCl	Bromazepam	Brompheniramine	Papaverine
Bupropion HCl	JWH-01965-Hydroxyhexyl metabolite	Rabeprazole sodium	Cyclobenzaprine	Doxepin	Cholesterol
Buspirone HCl	JWH-073 (Spice Cannabinoid)		Mephedrone	Nitrazepam	Aspartame
Butalbital			Pregabalin	Difunisal	
cocaine Hcl					

BIBLIOGRAPHY

1. GREENE, S.A. and THURMON, J.C. (1988), Xylazine - a review of its pharmacology and use in veterinary medicine. *Journal of Veterinary Pharmacology and Therapeutics*, 11: 295-313.
2. Reyes, J.C., Negrón, J.L., Colón, H.M. et al. The Emerging of Xylazine as a New Drug of Abuse and its Health Consequences among Drug Users in Puerto Rico. *J Urban Health* 89, 519-526 (2012).
3. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 2nd ed. Davis: Biomedical Publications; 1982.
4. Stewart DJ, Inaba T, Lucassen M, Kalow W. *Clin. Pharmacol. Ther.* April 1979; 25 ed: 464, 264-8.
5. Ambre J. *J. Anal. Toxicol.* 1985; 9:241.