

Disposable
Circular Stapler (Tilt Top)




EndoGuidance[®]
CIRCULAR STAPLER (Tilt Top)

YES, WE HAVE IT. IT'S IN OUR NAME[®]

Surgical Division

*More **Products**
More **Solutions**
Solid **Results***






EndoGuidance
CIRCULAR STAPLER (Tilt Top)


Disposable

Circular Stapler (Tilt Top)



Key Advantages

Reversible Nail Seat



Integrated Locking and Guiding Structure

Tilt Top Circular Stapler with independent intellectual property rights

Specification/Model

Specification Model	Height of staples (mm)	Outer Diameter of staples box head (mm)	Diameter of ring knife (mm)	Number of staples
RJCS-ST21	4.5	21	12.4	18
RJCS-ST23	4.8	24	14.4	18
RJCS-ST26	4.8	26	16.4	20
RJCS-ST29	4.8	29	19.4	24
RJCS-ST32	4.8	32	21.4	30
RJCS-ST34	5.0	34	23.4	32



(888) 863 -1112 Ext. 422

6831 Edgewater Commerce Pkwy Ste 1101
Orlando, FL 32810

MedicalDisposables.us
SurgicalDivision@MedicalDisposables.us

HAVE QUESTIONS? PLEASE CALL (888) 863-1112