Disposable Full-Detachable Trocar



## **Surgical Division**

More **Products**More **Solutions**Solid **Results** 





### Disposable

## **Full-Detachable Trocar**

#### Double Seal (Latex Free)

Ensuring biocompatibility and safety Accommodating instruments of various sizes ranging from 5mm - 15mm



#### **Primary Seal**

Funneled entry, keeping airtightness and smoothness while instruments use or exchange



#### **Secondary Seal**

Cross-entry guarantees airtightness when instruments are not in use



#### Full-detachable Design

To facilitate specimen removal and realizing fast discharge during the surgery



#### **Rotation Device**

180° rotation stopcock conveniently operated by both hands to control insufflation and deflation



#### **Ergonomic Handle**

Arc is more comfortable for operating streamlined design



#### **Detached Secondary Seal**

Rotary disassembly design realizes fast detachment for convenient surgeons during the procedure



#### **Barb Thread Design**


Providing stronger fixation to the abdominal wall



Disposable

# Auto-locking Full-detachable Trocar

Bladed, Safety Shield




#### **Auto-locking Trocar**

Item No.	Size
MDK5-70	5mm x 70mm
MDK5-100	5mm x 100mm
MDK5-150	5mm x 150mm
MDK10-100	10mm x 100mm
MDK10-150	10mm x 150mm
MDK11-100	11mm x 100mm
MDK11-150	11mm x 150mm
MDK12-100	12mm x 100mm
MDK12-150	12mm x 150mm
MDK15-100	15mm x 100mm
MDK15-150	15mm x 150mm

## **Bladeless Full-detachable Trocar**

Bladeless, Less Trauma




#### **Bladeless Trocar**

Item No.	Size
MD5-70	5mm x 70mm
MD5-100	5mm x 100mm
MD5-150	5mm x 150mm
MD10-100	10mm x 100mm
MD10-150	10mm x 150mm
MD11-100	11mm x 100mm
MD11-150	11mm x 150mm
MD12-100	12mm x 100mm
MD12-150	12mm x 150mm
MD15-100	15mm x 100mm
MD15-150	15mm x 150mm



### Disposable

# Optical Full-detachable Trocar Optical, Safety Shield



#### **Optical Trocar**

ltem No.	Size
MDO5-70	5mm x 70mm
MDO5-100	5mm x 100mm
MDO5-150	5mm x 150mm
MDO10-100	10mm x 100mm
MDO10-150	10mm x 150mm
MDO11-100	11mm x 100mm
MDO11-150	11mm x 150mm
MDO12-100	12mm x 100mm
MDO12-150	12mm x 150mm
MDO15-100	15mm x 100mm
MDO15-150	15mm x 150mm

## Hasson Full-detachable Trocar

Blunt-round Tip, Less trauma



-Hasson trocar is most suitable for obese/abdominal adhesion patients, which can greatly reduce the trauma area and herniation rate

#### **Hasson Trocar**

Item No.	Size
MDH10-10	10mm x 10mm
MDH12-10	12mm x 10mm



Like Us!







Follow Us!



& (888) 863 -1112 Ext. 422

© 6831 Edgewater Commerce Pkwy Ste 1101 Orlando, FL 32810

the penetration

MedicalDisposables.us

SurgicalDivision@MedicalDisposables.us