Disposable
Suction Irrigation System



**Surgical Division** 

More **Products**More **Solutions**Solid **Results** 





Disposable

# **Suction Irrigation System**

Integrated with Smoke Filter Type





Disposable

# **Suction Irrigation System**

Integrated with Smoke Filter Type

### **Ergonomic Handpiece**

Ambidextrous design allows it to be used with left or right hand comfortably

### **Completely Clear PC Material**

Provide the best view of the working channels of the handpiece part

#### **Color-coded Buttons**

For simple and quick distinction between suction and irrigation buttons

## **Intuitive Protrusion Design**

Convenient, suction button with protrusion design, no need to look to distinguish



## **Specimen Collection Portal**

Built-in sealing ring design ensures air tightness while enabling fast sampling

The integrated design of the lid and the sample collection channel is convenient for surgeons to open and close quickly, and is not easy to lose



Disposable

# **Suction Irrigation System**

Integrated with Smoke Filter Type



#### **Smoke Filter**

Remains pneumoperitoneum when the smoke filter is on Maintains visibility, improving overall efficiency in the operating room During surgeries reduces adverse effects on surgeons and patients.

### Specification/Model

Standard	With Smoke Filter	Size
MDSI-33	MDSI-33S	5 mm x 330 mm
MDSI-45	MDSI-45S	5 mm x 450 mm











