
Disposable Veress Needle



**Surgical Division** 

More **Products**More **Solutions**Solid **Results** 





# Laparoscopic

# Veress Needle

The cutting-edge Veress Needle revolutionizes minimally invasive surgeries, prioritizing precision and safety. Medical professionals can now perform procedures with unmatched accuracy and minimal patient discomfort.

#### **Enhanced Safety:**

Advanced safety features reduce the risk of accidental injury during insertion.

# **Optimal Precision:**

The finely tapered design enables precise entry into the body cavity, ensuring accuracy in targeted procedures.

### Versatility:

Caters to various minimally invasive procedures.

# Easy Handling:

Ergonomically designed to facilitate seamless maneuverability during surgery.

### Proven Efficiency:


Backed by extensive clinical research and positive outcomes, the Veress Needle has garnered trust among medical professionals worldwide.



Big size spring-loaded stylet provides loud audible "click" when passing through tissue layers



Visible red indicator indicate the needle position for safe insertion



Product Code	MDVN-120	MDVN-150
Size	120mm length	150mm length
25 Box		













