

Disposable
Electrosurgical Grounding Pad



EndoGuidance[®]
ELECTROSURGICAL PAD

YES, WE HAVE IT. IT'S IN OUR **NAME**[®]

Surgical Division

*More **Products**
More **Solutions**
Solid **Results***





EndoGuidance
ELECTROSURGICAL PAD

Disposable

Electrosurgical Grounding Pad Series

Secure the Circuit. Protect the Patient

MD EndoGuidance Electrosurgical Grounding Pad, Designed for Safety and Trusted by Surgeons.

The MD EndoGuidance Grounding Pads complete the electrical circuit during electrosurgery, allowing energy to safely return to the ESU. If a fault is detected at the patient return electrode, the system automatically deactivates the generator to prevent injury. Engineered to minimize burn risks and enhance comprehensive surgical safety.

Strong Adhesive Backing

Ensures the pad stays securely in place throughout the entire procedure.

Split Pad Technology

Monitors skin contact in real-time, helping prevent patient burns.

Latex-Free Construction

Safe for patients with latex allergies—gentle on all skin types.

Adult & Pediatric Sizes

Sized for safety—perfect fit for every patient, every time.


Universal Compatibility

Designed to work seamlessly with most electrosurgical units.

Peel-and-Apply Design

Quick to apply, easy to use—saves setup time in the OR.

MD-SY-MW2



Medical Disposables Corp.
YES, WE HAVE IT. IT'S IN OUR NAME

Like Us!



Follow Us!



Follow Us!



(888) 863-1112

6831 Edgewater Commerce Pkwy Ste 1101
Orlando, FL 32810

MedicalDisposables.us

SurgicalDivision@MedicalDisposables.us

HAVE QUESTIONS? PLEASE CALL (888) 863-8112